

Clear-Bite Biopsy Forceps




Oval Type

○ Oval shaped with small hole

○ Classification of application area by color

Different 2 color handles for intuitive identification of use

- **Violet** : Upper GI, Bronchoscopy
- **Orange** : Lower GI



Specification

Model	Cup Outer Diameter (@)	Working Length	Sheath O.D
FM-EF0001	Φ 1.72mm	1,600mm	Φ 1.63mm
FM-EF0002		2,400mm	
FM-EF0003	Φ 2.3mm	1,600mm	Φ 2.02mm
FM-EF0004		2,400mm	
FM-EF0005	Φ 1.8mm	1,100mm	Φ 1.7mm

Packaging : 20ea / box

Clear-Bite Biopsy Forceps

Round Type




- **Round Shape with small hole**

Effective for bleeding control

- **Coated Sheath**

Improving straightness when inserting into the scope.



Specification

Model	Cup Outer Diameter (a)	Working Length	Sheath O.D
FM-EF0003(Coated)(Round)	Φ 2.3mm	1,600mm	Φ 2.3mm
FM-EF0004(Coated)(Round)		2,300mm	


Packaging : 20ea / box

Needle Type



- **Needle for anchorage**

- **Oval shaped with small hole**



Specification

Model	Cup Outer Diameter (a)	Working Length	Sheath O.D
FM-EFN004	Φ 2.3mm	2,300mm	Φ 2.02mm

Packaging : 20ea / box