


Endo-Bite Biopsy Forceps



○ Reusable biopsy forceps of oval type



Specification

Model	Cup Outer Diameter (a)	Working Length	Sheath O.D
FM-REF0001	Φ 1.72mm	1,600mm	Φ 1.63mm
FM-REF0002		2,400mm	
FM-REF0003	Φ 2.3mm	1,600mm	Φ 2.02mm
FM-REF0004		2,400mm	

Packaging : 5ea / box

Clear- CoaBite




○ Design for Hot Biopsy

- Active mode (Hot biopsy)
Small polyps can be removed with a high-frequency current
- Non-Active mode enable to obtain a clean tissue

○ Oval shaped cup with small hole

○ Insulated cup inside

Minimizes damage to the biopsy tissue.



Specification

Model	Cup Outer Diameter (a)	Sheath O.D	Working Length
FM-EB0001 A	Φ 1.72mm	1,600mm	Φ 2.35mm
FM-EB0002 A		2,400mm	
FM-EB0003 A	Φ 2.3mm	1,600mm	Φ 2.35mm
FM-EB0004 A		2,400mm	

Packaging : 5ea / box