

# Clear-Hemostat



### 360° Coagulation Head

The silver tip can maintain a constant temperature and prevent sticking to the mucous membrane.  
Effective hemostasis by pressing with coagulation head



### Irrigation Port

Convenient when using scope without water jet function  
Easy to identify bleeding site.

## Specification

Model	Probe Head O.D	Working Length	Feature
FM-EH0001	Φ 2.3mm	2,400mm	Disposable

Packaging : 5ea / box

# ClearCoajet



injection mode

coagulation mode



### Hybrid hemostatic device

Injector + coagulator type facilitates the procedure devices  
by reducing the time to change devices during the procedure

### Irrigation Port

Convenient when using scope without water jet function  
Easy to identify bleeding sites

### Back - Cut beveled design

Accurate and smooth injection is possible with the  
design of the two-stage processing (Back-Cut Bevel) at the  
end of the needle.

## Specification

Model	Type	Thickness	Length	Gauge	Working Length
FM-CJ00C2(1800)	Coil Type	0.1T	4mm	23G	1,800mm
FM-CJ00C2(2300)					2,300mm

Packaging : 5ea / box