


Clear-Retriever

Hexagonal Type



- Hexagonal - shaped loop**
 Loop design that maintains a stable shape even after removing foreign bodies or ingested objects several times
- Soft & Durable Mesh**
 Polyester material to reduce tearing and abrasion
 Soft mesh material minimizes tissue damage
- 1:1 rotation with the handle**


Specification



Model	Mesh size(a)	Sheath O.D	Type of Rotation	Working Length
FM-ER0008A	30mm	Φ 2.5mm	Rotation	2,400mm


Packaging : 10ea / box

Oval Type



- Oval - shaped loop**
 Useful to remove removed polyps or foreign bodies or ingested objects with standard loop design
- Soft & Durable Mesh**
 Polyester material to reduce tearing and abrasion
 Soft mesh material minimizes tissue damage
 • **Voluminous pouch** : Suitable for collecting foreign substances of various sizes with deep mesh depth
- 1:1 rotation with the handle**

Specification



Model	Mesh size(a)	Sheath O.D	Type of Rotation	Working Length
FM-ER0004A	30mm	Φ 2.5mm	Rotation	2,400mm

Packaging : 10ea / box