


ClearEndoclip



Ring band type fastening method

The ring band fastening way provides high fixing force for effective bleeding control.

·The thin rear part that helps secure the endoscopic view : The rear mount is thin, so it is advantageous to secure the endoscopic view even after fastening the clip

Product design for safe use

Disposable all-in-one clip with no risk of infection
 Reduced installation failure rate compared to existing reused products
 Minimize damage to the endoscope channel as the clip is protected by an outer tube.

1:1 rotation with the handle

Accurate bleeding site positioning and easy targeting

Specification



Model	Opening Width (a)	Angle	Working Length
FM-CLD11-090(1650)	11mm	90°	1,650mm
FM-CLD11-090(1850)			1,850mm
FM-CLD11-090(2300)			2,300mm

Packaging : 5ea / box