

# 02


## EMR & Polypectomy



ClearGrasp Snare




Clear-Retriever



## ClearGrasp Snare

### Oval Type



#### Oval - shaped loop

The most commonly used Oval type loop design


#### Various loop sizes and wire stiffness

· Loop size : 6mm~30mm

· Stiffness : Soft, Moderato, Stiff

#### Convenient handle design

· 2 types : The non-rotating and rotating handles can be selected according to the operator's convenience



STERILE




### Specification

Type	Model	Loop size (a)	Stiffness	Type of Rotation	Working Length
Oval Type	FM-ES0006(2400)(Stiff)A	6mm	Stiff	Rotation	2,400mm
	FM-ES0006(2400)(Soft)A		Soft		
	FM-ES0006(2400)(Moderato)A		Moderato		
	FM-ES0004(2400)(Stiff)A	10mm	Stiff		
	FM-ES0004(2400)(Soft)A		Soft		
	FM-ES0004(2400)(Moderato)A		Moderato		
	FM-ES0005(2400)(Stiff)A	15mm	Stiff		
	FM-ES0005(2400)(Soft)A		Soft		
	FM-ES0005(2400)(Moderato)A		Moderato		
	FM-ES0002(2400)(Stiff)A	24mm	Stiff		
	FM-ES0002(2400)(Soft)A		Soft		
	FM-ES0002(2400)(Moderato)A		Moderato		
	FM-ES0003(2400)(Stiff)A	30mm	Stiff		
	FM-ES0003(2400)(Soft)A		Soft		
	FM-ES0003(2400)(Moderato)A		Moderato		
	FM-ESN004(2400)(Stiff)A	10mm	Stiff	Non-Rotation	
	FM-ESN004(2400)(Soft)A		Soft		
	FM-ESN004(2400)(Moderato)A		Moderato		
	FM-ESN005(2400)(Stiff)A	15mm	Stiff		
	FM-ESN005(2400)(Soft)A		Soft		
	FM-ESN005(2400)(Moderato)A		Moderato		
	FM-ESN002(2400)(Stiff)A	24mm	Stiff		
	FM-ESN002(2400)(Soft)A		Soft		
	FM-ESN002(2400)(Moderato)A		Moderato		
	FM-ESN003(2400)(Stiff)A	30mm	Stiff		
	FM-ESN003(2400)(Soft)A		Soft		
	FM-ESN003(2400)(Moderato)A		Moderato		

Packaging : 10ea / box

# ClearGrasp Snare

## Crescent Type



### ○ Crescent - shaped loop


Easy to remove lesions in difficult locations  
Suitable for EMR, cap-assistant EMR, EMRC

### ○ Various loop sizes and wire stiffness

- Loop size : 10mm~24mm
- Stiffness : Soft, Stiff

### ○ Convenient handle design

- 2 types : The non-rotating and rotating handles can be selected according to the operator's convenience




## Specification


Type	Model	Loop Size(ⓐ)	Stiffness	Type of Rotation	Working Length
Crescent Type	FM-ES0004(Crescent)(2400)(Stiff)A	10mm	Soft	Non-Rotation	2,400mm
	FM-ES0004(Crescent)(2400)(Soft)A		Stiff		
	FM-ES0005(Crescent)(2400)(Stiff)A	15mm	Soft		
	FM-ES0005(Crescent)(2400)(Soft)A		Stiff		
	FM-ES0002(Crescent)(2400)(Stiff)A	24mm	Soft		
	FM-ES0002(Crescent)(2400)(Soft)A		Stiff		

Packaging : 10ea / box

## Wire Thickness



## Cold type (Shield)



### ○ Cold (Shield - shaped) loop

Angled shield-shaped design for precise approach over the lesion

### ○ Cold & Hot polypectomy available

- Multi Function Plug : Optimized to perform for both hot and cold snare polypectomy.



### ○ Various loop sizes and wire stiffness

Suitable for removing small polyps under 5mm  
and small polyps 6~9mm

- Loop size : 6mm, 10mm, 15mm
- Stiffness : Soft (0.24mm, 0.3mm)

### ○ Convenient handle design

- 2 Types : The non-rotating and rotating handles can be selected according to the operator's convenience



## Specification

Type	Model	Loop size(㊤)	Stiffness	Type of Rotation	Working Length
Cold Type	FM-ES0006(Cold)(2400)A	6mm	Soft (0.24)	Rotation	2,400mm
	FM-ES0006(Cold)(2400)(0.3)A		Soft (0.3)		
	FM-ES0004(Cold)(2400)A	10mm	Soft (0.24)		
	FM-ES0004(Cold)(2400)(0.3)A		Soft (0.3)		
	FM-ES0005(Cold)(2400)A	15mm	Soft (0.24)		
	FM-ES0005(Cold)(2400)(0.3)A		Soft (0.3)		
	FM-ESN006(Cold)(2400)A	6mm	Soft (0.24)	Non-Rotation	
	FM-ESN006(Cold)(2400)(0.3)A		Soft (0.3)		
	FM-ESN004(Cold)(2400)A	10mm	Soft (0.24)		
	FM-ESN004(Cold)(2400)(0.3)A		Soft (0.3)		
	FM-ESN005(Cold)(2400)A	15mm	Soft (0.24)		
	FM-ESN005(Cold)(2400)(0.3)A		Soft (0.3)		

Packaging : 10ea / box