

○ **2 models available for the user's convenience**

○ **4.4Fr Tapered Tip Design**

Can be used for stenosis and small papilla, easy intubation

○ **Rounded Tip Forming**

Minimize tissue damage upon contact by round forming the tip of the tube



○ **Pre - Curved distal [11 o'clock]**

Provides cutting position between 11 o'clock to 12 o'clock

○ **3 Lumens**

Guidewire, cutting wire, 3 lumens for contrast injection

Color Marker



Specification

Model	Distal Tip Size (Fr)	Tip Length(㊟)	Cutting Wire Length(㊟)	Insulation Tube(㊟)	Working Length
FT-5201	4.4Fr	5mm	20mm	Yes	2,200mm
FT-5251			25mm		
FT-5301			30mm		
FT-7201		7mm	20mm		
FT-7251			25mm		
FT-7301			30mm		
FT-5200		5mm	20mm	No	
FT-5250			25mm		
FT-5300			30mm		
FT-7200		7mm	20mm		
FT-7250			25mm		
FT-7300			30mm		

Packaging : 1ea / box