

# Clear-Accu Biopsy Forceps



## Premium biopsy forceps

### Biopsy forceps with ultra-precision processing technology applied

Sharp cutting performance minimizes bleeding and enables clean tissue acquisition

- Swing cup helps easy access to lesions in difficult locations
- 2-Wire structure provides uniform force to the end of the cup

### Oval shaped with small hole

Oval-shaped cup with a small hole for obtaining more tissue biopsy samples


### Coated Sheath

Improving straightness when inserting into the scope

### Easy cup size identification by color

Different 2 color handles for intuitive identification of use

- **Violet** : 1.8mm
- **Orange** : 2.3mm



## Specification

Model	Cup Outer Diameter (@)	Working Length	Sheath O.D
FM-AF0001	Φ 1.8mm	1,600mm	Φ 1.7mm
FM-AF0002		2,300mm	
FM-AF0003	Φ 2.3mm	1,600mm	Φ 2.3mm
FM-AF0004		2,300mm	
FM-AF0005	Φ 1.8mm	1,200mm	Φ 1.7mm

Packaging : 20ea / box